

# Legacy® IPC promotes early discharges, enhances bed management, reduces readmissions due to DVT and improves safety scores!

**Problem:** Nurses are reluctant to discharge patients early since SCDs are not provided to patients for use throughout home recovery.

**Solution:** The transformational Legacy® IPC promotes early discharges and an At-Home DVT Prevention Therapy Program, proven to reduce readmissions due to DVT by 70%.



Kendall SCD™ 700

The incidence of inpatient DVT within the VA Hospital System is approximately 1.5%. However, readmissions due to DVT within 30-90 days of discharge is 2.8%-13% depending on patient's risk assessment and type of surgery. Hospitals performing 3,000 in-patient surgical procedures experience approximately 84 readmissions due to DVT, with 45% developing Post Thrombotic Syndrome (PTS).

**Legacy's At-Home DVT Prevention Program provides continuity of care from day of discharge and throughout home recovery until returning to an active lifestyle.**



Silent in Operation



Weighs 5 oz.



Promotes mobility



High compliance



Wear at work

## Steps Involved in Legacy's Early Discharge and At-Home DVT Program:

1. Day of discharge, nurses remove SCDs and apply the Legacy® IPC.
2. Nurses educate patients on how to apply and wear the Legacy® IPC.
3. Nurses stress the importance of wearing Legacy® IPC for 18 hours daily to prevent DVT.
4. Patients leave hospital wearing Legacy® IPC which reinforces the importance of compliance.
5. Ease of use, wearability, and comfort drives patient compliance above 85%.
6. Patients prefer Legacy's At-Home DVT Prevention Program by 90% versus extended hospital stay.

### The American College of Chest Physicians (ACCP) recommend:

- ☐ 18 hours of daily DVT prevention therapy.
- ☐ Daily compliance of at least 75%.
- ☐ Device be portable and promotes wearability.
- ☐ Battery-powered IPC device.

**The Legacy® IPC is a Game Changer in Preventing Readmissions Due to DVT Developing During Home Recovery!**

## Providing Veterans with Transformational Technology Prevents DVT and a Lifetime in Suffering with Post Thrombotic Syndrome.

Early ordering ensures Veterans receive Legacy® IPC on the day of discharge.



Minimum storage space



Custom application



One size fit all

### Legacy® IPC, is the only portable device evaluated at a major U.S. hospital: Mayo Clinic

- ☐ Trial involved 280 patients for four-months.
- ☐ Daily patient compliance above 85%.
- ☐ Average patient wear time was more than 18 hours daily.
- ☐ 0% hospital-acquired DVT cases.
- ☐ No cases involving falls and/or bleeding.

### Why provide the Legacy® IPC device to Veterans?

- ☐ Legacy's Piezo pump weighs only 5 grams that promotes ambulation.
- ☐ Piezo pump generates the same inflation pressure as bulky SCDs.
- ☐ Studies show 45 mmHg of compression therapy every minute for 15 seconds prevents stasis.
- ☐ No cumbersome tubes or power cords causing accidental falls and injuries.
- ☐ Open-cell fabric of the Legacy Comfort Sleeves promote breathability and real comfort.
- ☐ Legacy is the only DVT device that operates in total silence, which promotes quality sleep.
- ☐ Legacy's Lithium battery provides 24 hours of uninterrupted therapy on a full charge.
- ☐ Transformational technology promotes high compliance and positive outcomes.

### VasoCARE's FSS Contract 36F79726D0055 (Jan. 15, 2026 – Jan. 14, 2031)

#### ☐ Legacy IPC Portable DVT Device System – \$477.00

Item # 001-000 includes:

- 2 Legacy Pumps
- 2 Bladders - Left and Right
- 2 Comfort Sleeves - Left and Right
- 1 Dual USB Lead
- 1 Power Supply

**Legacy® IPC Offers a Drug-Free Option in Preventing DVT and Saving Lives!**

**VasoCARE** LLC

6554 Florida Blvd, Suite 123, Baton Rouge, LA 70806

(O) 225-924-0444 (F) 866-455-5150 Email: [orders@vasocare.com](mailto:orders@vasocare.com) [www.vasocare.com](http://www.vasocare.com)