

VasoCARE^{LLC}

Tomorrow's Technology for
Today's Medical Needs

Overview

Legacy[®] IPC (Intermittent Pneumatic Compression) is the first wearable DVT prevention device that promotes patient compliance above 85% and prevents hospital re-admissions due to DVT developing while Veterans are recovering at home.

Legacy[®] IPC is proven to improve patient compliance by 85% and reduce DVT by 70%.

Legacy[®] IPC offers a new level of care in pre-venting DVT.

Legacy[®] inflates the calf sleeves for 15 seconds at 45 mmHg every minute before deflating the sleeves for 45 seconds, allowing for maximum venous refill.

Dynamically emptying and refilling of the deep leg veins is key in preventing stasis and DVT.

Legacy[®] IPC is perfect for Veterans undergoing joint replacement surgery.

CORE DATA:

NAICS: 339113, 621999
Surgical Appliances



COMPANY SNAPSHOT:

CAGE Code: 75U54
Unique Entity: **FAEYHALH6FZ6**
FSS Contract#: **36F79726D0055**
Point of Contact: John Berges
Email: john@vasocare.com
Phone: 800-256-9979
Fax: 866-455-5150
Address: 6554 Florida Blvd.
Suite 123
Baton Rouge, LA 70806
Accreditation: HQAA Networks:
Medicare, Medicaid, BCBS,
Humana, UHC, etc.

KEY FEATURES:

- One button operation
- Compact/fully portable
- No cumbersome tubes/cords
- Sealed rechargeable battery
- Promotes mobility
- Breathable, open cell fabric
- 3-point system prevents slipping

TECHNICAL SPECIFICATIONS:

- Size: Pump is 4.1" x 1.7" x 0.67"
- Weights: 70 grams
- Mode of Operation: Cyclic
- Power: Lithium battery pack
- 24 hours of uninterrupted therapy

The Legacy[®] IPC System Produces Positive Outcomes!




Generational technology. One-size reduces inventory. Promotes ambulation.

Legacy[®] pump weighs 70 grams, yet mobilizes venous volume and velocity equal to big, bulky SCD units!

Legacy® IPC Wearable Device Technology

Legacy® IPC provides VTE prevention during hospital stay and throughout home recovery. No need for nurses removing and reapplying SCD limb sleeves throughout the day. Nurses find legacy® IPC superior to SCDs by over 85%. Legacy IPC promotes early discharge, prevents HA-DVT, and proven to prevent hospital readmissions due to DVT developing while Veterans are recovering at home.

Product	HCPC#	Description	Qty.	Price
Legacy® IPC	E0676	DVT prevention device	1	\$477.00
VascuEase System contains:				
2 Pumps - 2 Comfort Sleeves - 2 Bladders - 1 Power Supply - 1 Dual USB Lead - Mobile APP				
Contact Information:				
VasoCARE, LLC, 6554 Florida Blvd, Suite 123, Baton Rouge, LA 70806				
Office: 800-256-9979 or 225-924-0444 Fax: 866-455-5150				




Risk extends for 8 weeks.



3-Point System Prevents Slipping.



Advanced Pump Technology!



Legacy® Mobile APP.



Legacy® maximizes storage. Reusable system reduces costs.



Reduces post-op edema and accelerates healing!

Legacy® IPC is perfect for Veterans undergoing joint replacement surgery to wear throughout hospital stay and during post-surgical recovery while at home until fully ambulatory.

VasoCARE^{LLC}
 Tomorrow's Technology for
 Today's Medical Needs



VasoCARE, LLC
 is proud to serve
 our Veterans