

It's about protecting our Veterans from developing DVT while recovering at home. It's about saving lives and significantly reducing costs!

CAPABILITY STATEMENT www.vasocare.com

Problem: Nurses are reluctant to discharge patients early since SCDs are not provided to patients for use throughout home recovery.

Solution: The transformational Legacy[®] IPC promotes early discharges, provides high level DVT prevention therapy for Veterans recovering at home, and proven to reduce hospital readmissions due to DVT by 70%.

Legacy[®] IPC is a single-patient device that is used during home recovery until returning to an active lifestyle.

The ultra-portable device with Comfort Sleeve weighs only 5 ounces and powered by a Lithium battery that provides 24 hours of uninterrupted therapy on a single charge.

Legacy's Piezo pump weighs only 5 grams, delivers 45 mmHg for 12-15 sec. every minute, which is equal to bulky SCDs.

The open-cell fabric of the Comfort Sleeves promote breathability and real comfort.

Legacy is the only DVT device that operates in total silence.

It's about promoting high compliance and positive outcomes.

Unique Entity ID

FAEYHALH6FZ6

Legacy[®] IPC significantly reduces hospital costs, the incidence of re-admissions due to DVT by 70%, and prevents the development of Post Thrombotic Syndrome.

Legacy[®] IPC is the only portable device evaluated at major hospital: Mayo Clinic

- Trial involved 280 patients for 4 months.
- Daily patient compliance above 85%.
- Average wear time over 18 hours daily.
- 0% hospital-acquired DVT cases.
- No cases involving falls and/or bleeding
- Patients prefer Legacy's At-Home DVT Prevention Program by 90% versus an extended hospital stay.
- Legacy[®] IPC promotes mobility.

VasoCARE's FSS Contract 36F79726D0055 (January 15, 2026-January 14, 2031)

Legacy IPC Portable DVT System - \$477.00

Item # 001-000 includes:

- 2 Legacy Pumps
- 2 Bladders – Left and Right
- 2 Comfort Sleeves – Left and Right
- 1 Dual USB Lead
- 1 Power Supply

COMPANY SNAPSHOT

CAGE Code: 75U54
339113 Surgical Appliance & Supplies
621999 Miscellaneous Ambulatory Health Services

Government Business Point of Contact

John Berges

Email: john@vasocare.com

Phone: 225-924-0444

Fax: 866-455-5150

Address: 6554 Florida Blvd., Suite 123
Baton Rouge, LA 70806

FDA approval for the prevention of Deep Vein Thrombosis.

